



Vaccine Temperature Best Practices for **Frozen Vaccines–Fahrenheit (F)**

1 Store vaccine at ideal temperature



2 Record daily temperatures



Three Steps - Twice a Day: Temperatures should be checked and recorded first thing in the morning and before leaving at night.

- 1 Current Temperature:** The temperature that the freezer is right now.
- 2 Min/Max:** The coldest and warmest the freezer has been since you last reset the thermometer.
- 3 Reset:** The button you push after you have checked the Min/Max.



Best Practices


3 Take action if out of range!

- Contact your state or local health department immediately. Or if private vaccine call the manufacturer directly.
- Tell them the total amount of time the freezer was out of range.
- **Take your time** - Read and record temperatures accurately.
- **Make your mark!** Initial the log when recording temperatures.
- **Leave it blank** - if a temp was not recorded, leave the space blank!

Test Your Knowledge


Review the below temperature readings and select the correct answer.

1




- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

2




- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

3



- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

4



- A. Current temp and min/max are within range - no action necessary
- B. Current temp is within range, min/max out of range - take action
- C. Current temp is within range, min/max out of range - no action necessary
- D. Current temp and min/max are out of range - take action

5

Take action means (circle any that apply)

- A. Call the state/local VFC program (or manufacturer if private vaccine) for guidance
- B. Notify practice vaccine coordinator to change the temperature controls to get the freezer back in range
- C. Remove all vaccines that were out of range and discard them
- D. Discard any vaccine that does not look frozen